

Product Name :
Automatic Slide Staining Machine

Product Code :
NLAB-LABORATORYEM28005



Description :

Automatic Slide Staining Machine

Technical Specification :

Compact bench top model.

Perfect H & E Stainer for Routine Application.

High Specimen Through put.

Minimum Reagent Exposure.

Easy to clean.

Simple to operate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

