

Product Name :
Automatic Steel Knife Grinder

Product Code :
NLAB-LABORATORYEM28006



Description :

Automatic Steel Knife Grinder

Technical Specification :

This machine can grind all kinds of knife blade which length is below 240mm, and can replace the handcraft.

Technical Data:

Angle for adjustment of knife -edge: 20° -45° .

Interval time for turnover of knife: 5s, 10s, 25s

Grinding time: 0– 60 minutes set arbitrarily

Grinding knife size: Below 240 mm in length

Working voltage: AC220V 50Hz AC110V 60Hz

Power: 40 W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003