

**Product Name :**  
Automatic Sterilizer

**Product Code :**  
NLABMEDICAL-M32100020

**Description :**

Automatic Sterilizer India, Automatic Sterilizer manufacturer India, suppliers, Medical Automatic Sterilizer India, Hospital Automatic Sterilizer India

**Technical Specification :**

Temperature Sensor controls heating inside the chamber

Temperature Switch prevents from over-heating

Safety Pressure Switch maintains condition in complete safety

Error Messages for improper operation

Shut-Down Program for unsafe procedure

Auto self inspection

Double safe design of door

Two Types Available:

Vacuum Type: Avi-5510

Three pre-vacuum & one post-vacuum

Non-Vacuum Type Type: AVI-5511 (CE Approved)

Accessories:

Trays (30 x 17 x 2.2 cm)

4 pcs

Tray Rack

1 pcs

Drain hose with connector

1 pc

Clamp for adjusting two front legs

1 pc

Water filter

1 pc

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003