

Product Name :
Automatic TBA Device for Determining The Softening Point

Product Code :
NLAB-EDUCATIONTVET11004



Description :

Automatic TBA Device for Determining The Softening Point

Technical Specification :

Heating with magnetic stirrer according to EN 1427

230V / 50Hz / 1ph

Programmable test for water or glycerin bath

Large display 320x240 pixels, indicating in real time during the test: bath temperature, test status and agitator speed

Unlimited memory thanks to the SD card or external USB key port

Set of 2 rods with spherical ends for calibration verification

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com