

Product Name :
Automatic Tissue Processor with Microprocessor

Product Code :
NLAB-LABORATORYEM28008



Description :

Automatic Tissue Processor with Microprocessor

Technical Specification :

Feather Touch key pad

Digital display for adjusting time sequence

Battery backup to retain data

12 stage timing sequence can be programmed with each step of 1 minute to 16 hours 39 minutes (999 min.)

Delay action of 1 hour to 99 hours

2 Electronically Controlled wax bath 10°C to 80°C

1 or 2 Stainless Steel Basket

24 tissue containers (Cassettes)

1 or 2 Stainless steel Basket Rotors

10 without lip glass beakers (½ ltr, 1 ltr or 2 ltr.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003