

Product Name :
Automatic Tissue Processor

Product Code :
NLAB-EDUSCHOOLEQ29003



Description :

Automatic Tissue Processor

Technical Specification :

Automatic Tissue Processor

- Compact and versatile equipment for staining process atomization from stage one of fixation of tissue and dehydration process to paraffin wax saturation.
- Glass door with locking arrangement provides easy viewing and safety from tempering.
- Continuous agitation is provided by rotating S.S.
- Tissue basket to make process more effective.
- A superior thermostat maintains desired temperature of wax bath.

Accessories:

- 1 - S.S. Basket rotor
- Tissue basket
- **24** - S.S. Tissue capsules
- **10** - 1 liter glass beakers
- **11** - Beaker covers
- **2** - S.S. Wax baths each fitted with thermostat, **5-24** hrs
- Calibrated timing discs
- One-notch pliers for programming time discs
- It works suitably on **220V** ac **50Hz** single phase

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003