

Product Name :
Automatic Vicat Apparatus

Product Code :
NLAB-CEMETLAB130017



Description :

Automatic Vicat Apparatus

Technical Specification :

Conforming to EN 196-3, 13279-2, 480-2, ASTM C191, C187

Large size 4,3? touch screen color display

LAN port for direct connection to PC of a single unit or connection to a LAN hub for creating a network with up to 32 independent units all controlled by a single PC. 1 LAN cable is included

USB port for data storage on pen-drive (included)

Test procedures can be customized and stored to match user-defined requirements

Can incorporate an integrated graphic printer showing test result and setting time plot

Large test space with easy accessibility

Automatic calculation of initial and final setting time at programmable penetration depth limits

Wide range of accessories including EN and ASTM/AASHTO parts, in-water testing kit, needle cleaning device, integrated printer, probes for testing consistency and gypsum

Minimum penetrations rate: 10 seconds

Penetration measurement by encoder

Power: 50 W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003