

**Product Name :**  
Automobile Stethoscope Engineering Model and Training System

**Product Code :**  
NLAB-ENGINEERINGLB40005



**Description :**

Automobile Stethoscope Engineering Model and Training System

**Technical Specification :**

Automobile Stethoscope Engineering Model and Training System for engineering schools.

Automobile Stethoscope Engineering Model and Training System, technical teaching equipment for engineering - The mechanic's stethoscope makes it easy to identify piston slap, worn gears, faulty valves, water pump failure, damaged gaskets, defective bearings and other mysterious bumps and squeaks from the engine of the vehicle to the dashboard. Featuring a shockproof handle for comfortable operation, this mechanic's stethoscope reaches 11-3/4". Perfect for any shop or garage, this stethoscope is a tool any mechanic shouldn't be without.

- Identify noises in your car from the engine to the dashboard
- Two-piece probe reaches 11-3/4"
- Shockproof handle

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003