

Product Name :
Automotive Battery Load Tester

Product Code :
NLAB-MECHANICALAB20044



Description :

Automotive Battery Load Tester

Technical Specification :

For testing batteries, starter, alternator, voltage regulator, and related circuitry in 20 sec. test.

Heavy duty cable and copper plated clamps.

Stainless steel housing.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003