

Product Name :
Automotive Belt Tension Gauge

Product Code :
NLAB-MECHANICALAB20033



Description :

Automotive Belt Tension Gauge

Technical Specification :

Could be used to adjust the tension on the cam belts where the movement of the timing belt is measured by load.

Should be suitable for use on cam belts with various belt thickness.

Should be clear incremental marking with knurled knob for precise measurement.

Dual scale allows reading from top or bottom of tool.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003