

**Product Name :**  
Automotive Circuit Tester

**Product Code :**  
NLAB-MECHANICALAB20048



**Description :**

Automotive Circuit Tester

**Technical Specification :**

Power on type tests app. 6, 12, 24 V circuits

Length approx.135mm cable with insulated alligator clip.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003