

| | |
|---|---|
| Product Name : Automotive Cut section model of fully floating differential And rear wheel mechanism (working) | Product Code : NLAB-MECHANICALAB20068 |
|---|---|



Description :

Automotive Cut section model of fully floating differential And rear wheel mechanism (working)

Technical Specification :

This model is made out of original used heavy vehicle parts, suitably sectioned to show clearly the action of Differential gear box, such as pinion crown wheel in Differential and Brake drums in Rear axle, A Crank handle is provided to demonstrate the model. The model is mounted on a sturdy iron frame.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003