

**Product Name :**  
Automotive Cylinder Compression Tester For Petrol Engine

**Product Code :**  
NLAB-MECHANICALAB20045



**Description :**

Automotive Cylinder Compression Tester For Petrol Engine

**Technical Specification :**

Full 2 1/2" dial

Color-coded quadruple calibrations with app.0-300 psi, 21 kg/cm<sup>2</sup>, 21 bar and 2100 kpa

App. 6 heavy duty stem

App. 13 durable flexible hose

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003