

Product Name :
Automotive Diesel Engine Compression Tester

Product Code :
NLAB-MECHANICALAB20039



Description :

Automotive Diesel Engine Compression Tester

Technical Specification :

Test on both direct and indirect diesel injection engines.

Various size adaptors and quick release fitting for ease of use.

Easy to read large dual scale (psi/bar) gauge with 355mm flexible long hose aids easy use and quick release is fitted with pressure release valve for safe operation.

Packed in a rugged carry case.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003