

Product Name :
Automotive Diesel Injector Tester

Product Code :
NLAB-MECHANICALAB20040



Description :

Automotive Diesel Injector Tester

Technical Specification :

Should be used to adjust the injection nozzle opening pressure and to carry out leakage, spray pattern and chatter characteristics tests
For nozzle sizes R, S and T In basic design with 0 ... 40 MPa (0...400 bar) pressure gauge

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003