

**Product Name :**  
Automotive Digital Multimeter

**Product Code :**  
NLAB-MECHANICALAB20017



**Description :**

Automotive Digital Multimeter

**Technical Specification :**

**Specifications:**

- DC Voltage: 200m/2V $\hat{A}$  $\pm$ 0.5%, 20/200V/1000V $\hat{A}$  $\pm$ 0.8%
- AC Voltage: 200/750V $\hat{A}$  $\pm$ 1.2%
- DC Current: 200 $\hat{A}$  $\mu$ /2mA $\hat{A}$  $\pm$ 1.0%, 20m/200mA $\hat{A}$  $\pm$ 1.5%, 10A $\hat{A}$  $\pm$ 3.0%
- AC Current: 200/2K/200K/2M $\hat{A}$  $\pm$ 1.0%
- Diode Test: Forward DC Current: APPROX.1mA; Reversed
- DC Voltage: Approx.2.8V
- Transistor hFE: Test1-1000
- LCD display, easy readout; data hold.
- Max. Show Value: 1999 (3  $\hat{A}$  $\frac{1}{2}$  Display: 15mm LCD)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003