

**Product Name :**  
Automotive Gasoline Fuel Injection System Tester

**Product Code :**  
NLAB-MECHANICALAB20038



**Description :**

Automotive Gasoline Fuel Injection System Tester

**Technical Specification :**

Should be heavy-duty pressure tester kit includes a large dual scale (psi/bar) gauge, 4-ft hose, valve, 6-ft tubing with release valve, and various adapters Should be packed in a rugged case.

The adapters should be suitable for wide range of American, Japanese and European Vehicles.

The fuel pressure gauge could be hooked up or teed into system for ease of use.

**Specifications:-**

Temperature:-5~400C

humidity: = 95%

Air pressure: 70.0 ~106kPa

Power: AC 220V $\hat{\pm}$ 10%;50Hz $\hat{\pm}$ 1Hz

**Measurement scope:**

HC: 0 ~ 10000, 10<sup>-6</sup> (ppm) vol

CO: 0 ~ 10.0, 10<sup>-2</sup>(%)vol

CO<sub>2</sub>: 0 ~ 20.0, 10<sup>-2</sup> (%) vol

O<sub>2</sub>: 0 ~ 25.0, 10<sup>-2</sup>(%) vol

NO: 0 ~ 5000, 10<sup>-6</sup> (ppm) vol

oil temperature: 0 ~ 1200C

rotate speed: 0 ~ 10000rpm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003