

Product Name :
Automotive Grease Gun

Product Code :
NLAB-MECHANICALAB20030



Description :

Automotive Grease Gun

Technical Specification :

Capacity app.: 500 cc

Delivery pressure app.: 4500 psi / 310 bar, burst pressure app.: 10000 psi/ 690 bar

Finish: Chrome plated

Size app.: 352mmL x 62.5mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003