

Product Name :
Automotive Hand Held Tune-Up Analyzer

Product Code :
NLAB-MECHANICALAB20036

Description :

Automotive Hand Held Tune-Up Analyzer

Technical Specification :

Large display: 3-3/4 digit app. 4000 count, auto arranging meter with analog bar graph and additional setting symbol.

Record max. and min.

RS-232 interface to PC

Combine trig $\hat{A}\pm$ with duty cycle or pulse width

DC and AC Volts: 400mV, 4V, 40V, 400V, 500V.

Resistance: 400 Ω , 4K Ω , 40K Ω , 400K Ω , 4M Ω , 40M Ω .

Continuity buzzer

Diode test to check alternators, relays, solenoids.

K-type temperature probe should be included.

Currents DC and AC Current: 40mA, 400mA, 20A

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003