

Product Name :
Automotive Rack and Pinion Type

Product Code :
NLAB-MECHANICALAB20069



Description :

Automotive Rack and Pinion Type

Technical Specification :

This laboratory system is used to introduce combined drives and their correct assembly. The programme of exercises enables familiarisation with six different drive implementation and analysis methods: understanding the brief and the drawing, assembly, setting, adjusting, testing, and making calculations.

The flexibility of the set-up and the modularity of the components simplifies experimentation and implementation of the students' own ideas.

A robust tubular steel frame with a square profile and bearing elements provide the accuracy for the setting of precise gearing. All the system components are kept ready to hand and well protected in a housing system.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003