

Product Name :
Automotive Universal Radiator Pressure Tester

Product Code :
NLAB-MECHANICALAB20041



Description :

Automotive Universal Radiator Pressure Tester

Technical Specification :

Should be heavy duty pressure tester kit supplied with a heavy duty pump gauge unit, thermometer and multi adaptors.

Should be packed in a rugged case.

Heavy duty pump gauge unit features smooth action and easy to read color coded scale.

Should be multi-color coded aluminum radiator adaptors and quick release fitting with pressure release valve for ease of use

Suitable for wide range of vehicles applications: FIAT, Audi, BMW, Ford, Honda, Isuzu, Jaguar, Mazda, Mitsubishi, Suzuki, Toyota, Volvo etc.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

