

Product Name :
Axial Shoulder Exerciser

Product Code :
NLBMEDICAL-M3090001



Description :

Axial Shoulder Exerciser

Technical Specification :

Combo unit for shoulder & Supination – Pronation exercises, unit is used in the same way as a conventional shoulder wheel and rotary wrist machine. Unit is suitable to improve circulation, muscle strength, joint range of motion and coordination.

Technical :

Construction : Precision machined metal drum is fitted with calibrated sensitive resistance mechanism.

Resistance : Variable Calibrated Resistance, controllable from Zero to maximum by tightening the resistance strap.

Adjustment Bars : Two, made of square Stainless Steel tube with bright polish finish.

: One for Shoulder Exercise.

: One for Supination-Pronation movement.

Arc of Motion : Arc of Motion can be varied from 30cm to 90cm by fixing handle at the required place to record the improvement of shoulder movement.

Feedback : A 360 degree scale is provided on the drum to measure & record the degree of revolution from either direction (for Left or Right Shoulder) for immediate feedback to the patient.

Movement : Bi-directional operation for shoulder & wrist.

Height Adjustment : Drum Height is adjustable on 50cm long pipe for adjustment according to patient's need.

Mounting : Wall mountable unit, set include mounting hardware.

Finish : Powder Coated Finish.

We are Exporters and Suppliers of Axial Shoulder Exerciser at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003