

**Product Name :**  
Baby Bassinet Luxury Blue

**Product Code :**  
NLABMEDICAL-M2490005



**Description :**

Baby Bassinet Luxury Blue

**Technical Specification :**

Specification :-

Baby Bassinet Luxury is made of ABS Plastic.

The body section is mounted on 10 cms lockable castors.

It comes with the feature of Trendelenburg and reverse Trendelenburg positions.

It has plastic bassinet for clear visibility.

The Incline position is adjustable by Gas Spring upto 30 cms.

The height is adjustable by Gas spring upto 105 cms.

We are Manufacturer, Exporters and Suppliers of Baby Bassinet Luxury Blue at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003