

**Product Name :**  
Baby Bassinet Standard

**Product Code :**  
NLABMEDICAL-M2270001



**Description :**

Baby Bassinet Standard

**Technical Specification :**

Baby Bassinet Standard

Baby Bassinet is manufactured with heavy round pipes.

Consumables are kept at the Bottom Shelf.

Acrylic transparent Crib with mattress.

The body section is mounted on 5cm castors.

The crib has positions like Trendelenburg / Rev. Trendelenburg.

Product Size: 70 L x 45 W x 80 H cms.

Knock Down Construction to save freight.

We are Manufacturers, Exporters, and Suppliers of Baby Bassinet Standard at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003