

Product Name :
Baby Bassinet Standard

Product Code :
NLABMEDICAL-M2490003



Description :

Baby Bassinet Standard

Technical Specification :

Specification :-

Baby Bassinet is manufactured with heavy round pipes.

Consumables are kept at the bottom Shelf.

Acrylic transparent Crib with mattress.

Body section is mounted on 5cm castors.

Crib has positions like Trendelenburg / Rev. Trendelenburg.

Product Size: 70 L x 45 W x 80 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Baby Bassinet Standard at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003