

**Product Name :**  
Baby Bassinet with Hydraulic

**Product Code :**  
NLABMEDICAL-M2810002



**Description :**

Baby Bassinet with Hydraulic

**Technical Specification :**

Baby Bassinet with Hydraulic

Baby Bassinet is manufactured with heavy round pipes.

Consumables are kept at the Bottom Shelf.

Acrylic transparent Crib with mattress.

Panels are Foldable.

The body section is mounted on 5cm castors.

The crib has positions like Trendelenburg / Rev. Trendelenburg.

Product Size: 70 L x 45 W x 80 H cm.

Knock Down Construction to save freight.

We are Manufacturer, Exporters, and Suppliers of Baby Bassinet with Hydraulic at tender and export competitive prices from India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003