

Product Name :
Baby Bassinet with Hydraulic

Product Code :
NLABMEDICAL-M2490002



Description :

Baby Bassinet with Hydraulic

Technical Specification :

Specification :-

Baby Bassinet is manufactured with heavy round pipes.

Consumables are kept at the bottom Shelf.

Acrylic transparent Crib with mattress.

Panels are Foldable.

Body section is mounted on 5cm castors.

Crib has positions like Trendelenburg / Rev. Trendelenburg.

Product Size: 70 L x 45 W x 80 H cm.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Baby Bassinet with Hydraulic at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003