

**Product Name :**  
Baby Bassinet with Plastic Crib

**Product Code :**  
NLABMEDICAL-M2490001



**Description :**

Baby Bassinet with Plastic Crib

**Technical Specification :**

Baby Bassinet with Plastic Crib is manufactured with heavy round pipes.

Removable molded Plastic Crib without Joints.

The body section is mounted on 10cm castors, 2 with Brakes.

Trendelenburg Position Inclination up to 12 degrees.

Stainless Steel shelf at the bottom.

Height adjustable approx. 85 – 100 cm.

Product Size : 90 L x 50 W x 85-100 H cm.

We are Manufacturers, Exporters, and Suppliers of Baby Bassinet with Plastic Crib at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003