

Product Name :
Baby Bassinet with Utility Box

Product Code :
NLABMEDICAL-M2490004



Description :

Baby Bassinet with Utility Box

Technical Specification :

Specification :-

Frame is made of Round Tube.

The body section is mounted on 5cm castors.

It comes with IV Rod.

Acrylic transparent Crib with mattress.

The Crib has positions like Trendelenburg / Rev. Trendelenburg.

Wooden utility box below the Crib.

Product Size: 70 L x 45 W x 80 H cms.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Baby Bassinet with Utility Box at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003