

Product Name :
Bacteriological Incubator

Product Code :
NLABMEDICAL-M3170002

Description :

Bacteriological Incubator India, Bacteriological Incubator manufacturer India, suppliers, Medical Bacteriological Incubator India, Hospital Bacteriological Incubator India

Technical Specification :

Inner chamber made of aluminum or stainless steel sheet and outer made of mild steel sheet finished in powder coated color combinations.

Gap filled with glass wool for temp. Insulation.

Double walled door with double glass window for observation in the chamber.

Temperature range 5° C above ambient to 800C ± 10C controlled by a thermostat.

L-shaped thermometer provided along with pilot lamps, on-off switch temp.

Control knob etc.

It works on 220V AC 50Hz single phase.

Optional:

Air circulating fan

Timer 24 hours

Digital temp. Controller-cum-indicator

Digital temp. Indicator-cum-controller with computer data logger.

Chamber Size in mm

300 x 300 x 300

355 x 355 x 355

455 x 455 x 455

455 x 605 x 455

605 x 605 x 455

605 x 605 x 605

605 x 910 x 605

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003