

Product Name :
Ball Resilience Tester Machines

Product Code :
NLAB-TESTINGLABM92001



Description :

Ball Resilience Tester Machines

Technical Specification :

Technical Specification:

Diameter of steel sphere: 16 mm

Weight of steel sphere: 16.3 g

Dropping height: 500 mm

Height of test specimen: 50 to 150 mm

Inner diameter of acrylic tube: 40 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

