

Product Name :
Bandsaw Brazing Machine

Product Code :
NLAB-MACHINELABWK30017



Description :

Bandsaw Brazing Machine

Technical Specification :

Bandsaw Brazing Machine .

Indian

make Bandsaw Brazing Machine used for brazing the ends of bandsaw blades from 3 mm to 30 mm width capacity 3KVA suitable for operation on AC 220 volts single phase, 50 cycles electric supply with switch.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003