

Product Name :
Bar Breaking Apparatus

Product Code :
NLAB-EDUSCHOOLEQ152001



Description :

Bar Breaking Apparatus

Technical Specification :

The apparatus is supplied with a pack of ten breaking bars.

To show the forces which can be exerted during thermal expansion or contraction.

Comprising a heavy cast iron frame 340 x 120 x 100 mm with slotted end pillars to carry a stout iron bar.

The bar is threaded at one end for a large tensioning nut and has holes at the other end to accommodate the cast iron breaking bars.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003