

Product Name :
Bar Magnet

Product Code :
NLAB-ENGINEERINGLB12004



Description :

Bar Magnet

Technical Specification :

Bar Magnet Lab Equipments Manufacturer, Suppliers & Exporters

Made of steel, 12 x 5 mm in pairs ; half red & half blue (for N & S) spray painted, with keepers..

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003