


Product Name :
Barometer Aneroid Wall Mounting

Product Code :
NLAB-EDUSCHOOLEQ139001



Description :

Barometer Aneroid Wall Mounting

Technical Specification :

- Wall mounting.
- Dial 100 mm diameter, mounted in a plastic case 150 mm diameter.
- Graduated in inches and millibars

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003