


**Product Name :**  
Barometer Aneroid Wall Mounting

**Product Code :**  
NLAB-EDUSCHOOLEQ139001



**Description :**

Barometer Aneroid Wall Mounting

**Technical Specification :**

- Wall mounting.
- Dial 100 mm diameter, mounted in a plastic case 150 mm diameter.
- Graduated in inches and millibars

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003