

Product Name :
Barometer Demonstration Aneroid

Product Code :
LiquidsPressure0002



Description :

Barometer Demonstration Aneroid

Technical Specification :

- Encased in a transparent acrylic case which allows full view of the movement.
- The action can be demonstrated with the help of a rubber bellow which increases the pressure when squeezed.
- Supplied complete with rubber bellow.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773