


**Product Name :**  
Basic Metacentric Height

**Product Code :**  
NLAB-MECHANICALAB230002



**Description :**

Basic Metacentric Height

**Technical Specification :**

The advanced meta-center can be varied to produce stable and unstable equilibrium. The equipment consists of a plastic rectangular floating pontoon, the center of gravity of which can be varied by an adjustable weight, which slides and can be clamped in any position on a vertical mast. A single plumb bob is suspended from the mast, which indicates the angle of heel on a calibrated scale. A weight with lateral adjustment enables the degree of heel to be varied and hence the stability of the pontoon determined. The equipment does not require a separate water tank as it may be used

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003