

Product Name :
Basic PC Network System Trainer

Product Code :
NLAB-ENGINEERINGLB54052



Description :

Basic PC Network System Trainer

Technical Specification :

The components are fully compatible with most LAN software including Novell, Windows NT, Lotus Notes Windows 2000 or Linux etc.

The trainer comes with standard components and peripherals for a Network and provides all necessary modules and parts for building Thin net system (Thin Coax) and Twisted Pair Network (10/100-Base T).

PC Network System Trainer is designed for flexibility and simplicity in building LAN Systems.

We are well-known manufacturers, OEM suppliers of Basic PC Network System Trainer for Applied Machines. Contact us for high quality Basic PC Network System Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003