

Product Name :
Basket Stretcher

Product Code :
NLABMEDICAL-M2110004



Description :

Basket Stretcher

Technical Specification :

Specification:-

- * This stretcher is made with an eye on a very range of emergency situations, and it's sturdy and flexible for use in all types of aid in special situations.
- * Its reliable fittings enable first-aid personnel to operate rapidly and safely.
- * With its special sling equipment, the stretcher is ideal for lifting and transporting by helicopter.
- * It has an adjustable feet-securing mechanism, safety belts, and a mattress.
- * It is also designed to be secured inside the helicopter.
- * All the materials used have dimensions suitable to work and they are fireproof.
- * They don't release any toxic or polluting substances and they are protected to obtain high resistance to wear and corrosion.

Technical Parameter :

* Product Size (L x W x H) :- 216 x 61 x 19cm

* Packing Size (1 pc) :- 220 x 65 x 22cm

* Load Bearing:- 159 kg

* N.W (1 pc) :- 18 kg

* GW :- 20 kg

We are Exporters and Suppliers of Basket Stretcher at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003