

**Product Name :**  
Batch Distillation Pilot Plant

**Product Code :**  
NLAB-ENGINEERINGLB36018



**Description :**

Batch Distillation Pilot Plant

**Technical Specification :**

Tank for the collection of distillate with 1 litre capacity.

The tube in tube heat exchangers manufactured from stainless steel

6 digital thermometers with probes protected by stainless steel shield sleeves

Thyristor unit for controlling the power to the reboiler electric heater

Variable-area flow meter made of glass and stainless steel with a range 20 to 200 litres/hour

Emergency stop button

Distillation column height 600 mm with 7 sieve trays.

Reboiler with 3 litres capacity complete with 1600 W electric heater

Condenser with reflux head

2 programmable timers for the control of reflux ratio.

We are well-known manufacturers, OEM suppliers of Batch Distillation Pilot Plant for Chemical Engineering Lab Equipments. Contact us for high quality Batch Distillation Pilot Plant for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003