

**Product Name :**  
Beadbug Microtube Homogenizer

**Product Code :**  
NLAB-EDUSCHOOLEQ45043



**Description :**

Beadbug Microtube Homogenizer

**Technical Specification :**

Bead Bug is a high energy benchtop homogenizer that sets the performance standard for personal sized cell disruption and lysis instruments.

The BeadBug is the ideal homogenizing instrument for low to medium throughput laboratories.

The optimized mixing motion causes rapid cell disruption through constant high-velocity impact from the hardened micro-beads chosen specifically for your sample type.

Simultaneous homogenization of up to 3 samples takes place (often within 45 seconds) inside the disposable 2ml screw cap microtubes.

Developed for samples too tough to be homogenized in a common lab vortexer.

The BeadBug provides a convenient and economical alternative to the traditional, over-priced homogenizers on the market today.

With up to three tube capacity and a footprint of just 7 x 8".

Speed Range 3700 to 4000 rpm	
Capacity (Screw Cap / Non-Skirted)	
Timer: to 3 min.	
Dimensions 5.3 in./17 x 21 x 13.5 cm	
Weight 2.2kg	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003