

**Product Name :**  
Beakers Euro Design

**Product Code :**  
NLAB-EDUSCHOOLEQ122002



**Description :**

Beakers Euro Design

**Technical Specification :**

- Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance.
- These autoclavable Beakers have prominently raised graduations on the outer surface.
- Pouring a Solution becomes very convenient because of its tapered spout.

**Beakers**

25 Pcs

50 Pcs

50 Pcs

20 Pcs

50 Pcs

200 ml

2000ml

2000ml

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003