

Product Name :
Beam Apparatus

Product Code :
TN325



Description :

Beam Apparatus

Technical Specification :

Beam Apparatus

One of the ends can be fixed or given a known slope by applying a known moment at the end with the help of suspended loads. At the other end also a known moment can be applied. Vertical loads can be applied at various points along the span of the beam.

A dial gauge with 25mm travel (with a magnetic base) supplied with the apparatus. Apparatus to be supplied is complete with a supporting stand and a set of weights.

Spring Balance: 0 to 10N capacity

Weights: 10N/5N/2N

Beam: 19mm wide; 1300mm long

Cantilever support: Rigid pillar clamping arrangement

Dial Gauge: 0-25 mm; LC: 0.1 div

Bench Area: 2 x 06 Mts

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

