

Product Name :
Beam Moulds

Product Code :
NLAB-ENGINEERINGLB89002



Description :

Beam Moulds

Technical Specification :

Five sizes of beam molds are offered for casting concrete specimens for flexural strength testing.

These molds are available in cast iron & mild steel and are supplied complete with a base plate.

Optional Accessory:-

Tamping Bar, Mild Steel, 25mm x25mm square ramming face, 400mm long, 2Kg in weight.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003