

<b>Product Name :</b> Bed Full Fowler	<b>Product Code :</b> NLABMEDICAL-M2790006
<b>Description :</b> Bed Full Fowler India, Bed Full Fowler manufacturer India, suppliers, Medical Bed Full Fowler India, Hospital Bed Full Fowler India	
<b>Technical Specification :</b> Hospital Fowler Bed made of high quality materials, components and accessories  Four section sheet metal top  Fowler Bed with adjustable back section and knee-rest  Manually operated crank system for various positions  Made of high quality materials, components and accessories  Four section sheet metal top  Fowler Bed with adjustable back section and knee-rest  Manually operated crank system for various positions  Stainless steel handles for backrest adjustment  These fowler beds are equipped with easy lifting or collapsible side rails to ensure patient safety  Stainless Steel or ABS is used for the head and foot panel  Four section top of uniformly perforated CRCA sheet.  Backrest and knee rest section maneuvered by two crank mechanisms from the foot end.  ABS bows easily removable for access to the patient  Trolley type base mounted on 125 mm die  Size:	

Standard dimension:{LxWxH}

2030 x 900 x 600

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003