

**Product Name :**  
Bed Full Fowler

**Product Code :**  
NLABMEDICAL-M2790006

**Description :**

Bed Full Fowler India, Bed Full Fowler manufacturer India, suppliers, Medical Bed Full Fowler India, Hospital Bed Full Fowler India

**Technical Specification :**

Hospital Fowler Bed made of high quality materials, components and accessories

Four section sheet metal top

Fowler Bed with adjustable back section and knee-rest

Manually operated crank system for various positions

Made of high quality materials, components and accessories

Four section sheet metal top

Fowler Bed with adjustable back section and knee-rest

Manually operated crank system for various positions

Stainless steel handles for backrest adjustment

These fowler beds are equipped with easy lifting or collapsible side rails to ensure patient safety

Stainless Steel or ABS is used for the head and foot panel

Four section top of uniformly perforated CRCA sheet.

Backrest and knee rest section maneuvered by two crank mechanisms from the foot end.

ABS bows easily removable for access to the patient

Trolley type base mounted on 125 mm die

Size:

Standard dimension:{LxWxH}

2030 x 900 x 600

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003