

Product Name : Bed Side Locker-C	Product Code : NLABMEDICAL-M2730004
Description :	
Bed Side Locker-C	
Technical Specification :	
<p>The locker box with deep drawer is made of commercial board flushed with laminate, assembled on square chrome plated tubes</p> <p>The locker machine pressed from 20 gauge CRCA sheet</p> <p>All three sides closed having one box and a drawer</p> <p>Stainless Steel top</p> <p>Mounted on 5 cm DIA. Swivelling castors</p> <p>Finish in epoxy powder coated.</p> <p>Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities</p> <p>Bedside lockers optionally available without draw under the top</p> <p>A high quality material is used for the Bedside Cabinets</p> <p>These Bedside Lockers are helpful in keeping the medicines and other useful things in an organized manner</p> <p>Moreover, these products are ideally placed at the bedside so that all the stored products remain within the reach</p> <p>Side locker made with heavy duty CRCA MS sheet, duly machine pressed.</p> <p>Size Top :(40 cm sq. x 82 H) cm</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003