

Product Name : Bed Side Locker-D	Product Code : NLABMEDICAL-M2730005
Description :	
Bed Side Locker-D	
Technical Specification :	
Enamelled sheet steel frame	
Full length body	
One piece drawer	
Laminated drawer face and top	
Interior fitted with on steel shelf and perforated bottom	
The locker box with deep drawer is made of commercial board flushed with laminate, assembled on square chrome plated tubes	
The locker machine pressed from 20 gauge CRCA sheet	
All three sides closed having one box and a drawer	
Stainless Steel top	
Mounted on 5 cm DIA. Swivelling castors	
Finish in epoxy powder coated.	
Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities	
Size:- (40 cm sq. x 82 H) cm.	

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments