

Product Name :
Bedside Locker ABS

Product Code :
NLABMEDICAL-M2500006



Description :

Bedside Locker ABS

Technical Specification :

Specification :-

Bedside Locker is made of ABS Plastic.

It comes with two towel rail.

The body section is mounted on 5cm castors.

It comes with cupboard box.

It comes with one drawer for storage.

Product Size : 45 L x 45 W x 75 H cms.

We are Manufacturer, Exporters and Suppliers of Bedside Locker ABS at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003